

SEQUENCE LISTING

<110> Sinha, Sudhir K
Walker, Jerilyn A
Batzner, Mark A

<120> Assay for Quantitation of Human DNA Using Alu Elements

<130> P56885

<160> 9

<170> PatentIn version 3.1

<210> 1

<211> 19

<212> DNA

<213> Artificial

<220>

<223> Primer (Alu3) for inter-Alu PCR

<400> 1

gatcgcgcca ctgcactcc

19

<210> 2

<211> 20

<212> DNA

<213> Artificial

<220>

<223> Primer (Alu5) for inter-Alu PCR

<400> 2

ggattacagg cgtgagccac

20

<210> 3

<211> 22

<212> DNA

<213> Artificial

<220>

<223> The forward intra Yb8 primer for intra-Alu PCR

<400> 3

cgaggcgggt ggatcatgag gt

22

<210> 4
 <211> 20
 <212> DNA
 <213> Artificial

 <220>
 <223> The reverse intra Yb8 primer for intra-Alu PCR

 <400> 4
 tctgtcgccc aggccggact 20

 <210> 5
 <211> 21
 <212> DNA
 <213> Artificial

 <220>
 <223> The forward intra Yd6 primer for intra-Alu PCR

 <400> 5
 gagatcgaga ccacggtgaa a 21

 <210> 6
 <211> 19
 <212> DNA
 <213> Artificial

 <220>
 <223> The reverse intra Yd6 primer for intra-Alu PCR

 <400> 6
 tttgagacgg agtctcgtt 19

 <210> 7
 <211> 287
 <212> DNA
 <213> Homo sapiens

 <400> 7
 cgccggggcgc ggtggcgcac gcctgtaatc ccagcacttt gggaggccga ggcggggcgga 60
 tcacgaggtc aggagatcga gaccatcctg gctaacacgg tgaaaccccg tctctactaa 120
 aaatacaaaa aattagccgg gcgtggtggc gggcgccctgt agtcccagct actcgggagg 180

ctgaggcagg agaatggcgt gaacccggga ggcggagctt gcagtgagcc gagatcgcg	240
cactgcactc cagcctgggc gacagagcga gactccgtct caaaaaa	287

<210> 8
 <211> 294
 <212> DNA
 <213> Homo sapiens

<400> 8	
cgccggggcgc ggtggcgcac gcctgtaatc ccagcacttt gggaggccga ggcgggtgga	60
tcatgaggtc aggagatcga gaccatcctg gctaacaagg tgaaaccccg tctctactaa	120
aaatacaaaa aattagccgg gcgcggtggc gggcgctgt agtcccagct actcgggagg	180
ctgaggcagg agaatggcgt gaacccggga agcggagctt gcagtgagcc gagattcgcg	240
cactgcagtc cgcagtccgg cctggggcgac agagcgagac tccgtctcaa aaaa	294

<210> 9
 <211> 275
 <212> DNA
 <213> Homo sapiens

<400> 9	
cgccggggcgc ggtggcgcac gcctgtaatc ccagcacttt gggaggccga ggcgggcgga	60
tcacgaggtc aggagatcga gaccacggtg aaaccccgtc tctactaaaa atacaaaaaa	120
ttagccgggc gcagtggcgg gcgcctgtag tcccagctac tcgggaggct gaggcaggag	180
aatggcgtga acccggaagg cggagcttgc agtgagcgga gatcgcgcca cagcactccc	240
gcctgggcga cagaacgaga ctccgtctca aaaaa	275